



Recombinant *Drosophila simulans* Adenosine monophosphate-protein transferase FICD homolog (GD23409)

Product Code	CSB-CF455796DMJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B4Q4M7
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila simulans</i> (Fruit fly)
Sequence	<p>MCMEAEPSPPAQQQEQVNPPLCKAQNPKPARYRLVLLFVAGSLAAWTFHA LSSTNLVW KLRQLHHLPTAHYLQTRDEFALYSVEELNAFKEFYDKSVSDSVGASYTEAEQT NIKEALG ALRMAQDLYLAGKDDKAARLFEHALALAPRHPEVLLRYGEFLEHNQRNIVLAD QYYFQAL TISPSNSEALANRQRTADVQSLDERRLESLSKRDALSIAIHESNGALRRAKK EAYFQHI YHSVIGIEGNTMTLAQTRSILETRMAVDGKSIDEHNEILGMDLAMKYINASLVQKI DITIK DILELHRRVLGHVDPIEGGEFRRNQVYVGGHIPPGGDLALLMQRFERWLNSE HSSTLHP VNAAALAHYKLVHIHPFVDGNGRTSRLLMNTLLMRAGYPPV IIPKQQRSKYYHF LKLANE GDIRPFVRFIADCTEKTLDLYLWATSDLPQQIPMLIQTESEAGERLAQMSPNV AQRSSI LEFYESGSGDIP</p>
Source	in vitro E.coli expression system
Target Names	GD23409
Protein Names	Recommended name: Adenosine monophosphate-protein transferase FICD homolog EC= 2.7.7.n1
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.