



Recombinant Drosophila melanogaster Protein odr-4 homolog (CG10616)

Product Code	CSB-CF896148DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VTX8
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MRTVLLSKHDELYLEKCAQENQFSYGIIVGHQADLTKSVVVHLARNNEEDADL EDLSEVR LTISDINSQALASQWLSASKMCPGSFDVIGIFVSSVRSDDVNEQSAEFKNAKKL FSDIYD LLLKSNSSFGVYTTDIAQTDFVFLSYSLADKKVLCKNYSYGNNGTFTNMEFRFV DKPFEW IQLECSYDFDDVLPILDSSRRVNIEDQFQSMIVSVRKNLLASEVFLQNEVVETDI DLQAY IKKKKTKVDKLQPTSTTGGTATASSNTTDSLRLASEGIIGGTETIRASIVLPMKC QLSK PTDIKVREFSGTLHMSGIITSKVFCNPRNSIADVCRFLRDDVLRSLITRIQVYCD GLTDP YVTNEALYISEPPRRVFFSLPSEGPSASVGAVVQFSEYLFRGEAPTVVVAQAK QILDVDL DPETISVEAEGLPDDTHFNNCKMDADCIDDSGIMTSSMPKPELSRSLYMGVIA VALLVLL SSVALHFVLAER
Source	in vitro E.coli expression system
Target Names	CG10616
Protein Names	Recommended name: Protein odr-4 homolog
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.