



# Recombinant *Mentha piperita* Cytochrome P450 71A8 (CYP71A8)

<b>Product Code</b>	CSB-CF671473MQY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q42716
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Mentha piperita</i> (Peppermint) ( <i>Mentha aquatica</i> x <i>Mentha spicata</i> )
<b>Sequence</b>	MYSIAFMLVLRKMDEIISHTLAFQALVSLILLISITKWLSNSPKNKNSSPPSPRKL PILG NLLQLGSLPHHNLRSMARKHGPIMLLHLGSRPVSSRRRPRGNHENSRSRLR RPRGSRSA ALQLQGRVGGYGEYWRQLKTICVVQLLSNKRVSFRSVREEETELLMKKIGD SSGNVNLS HMFTQLTNDVVCRSAIGRKYGAGDENGKFLFLEILREFLELLGAISIGDFVPSLW WINRIN GFDRRVDRIAKEMDEFLEKVIHERLENPAAKAEENFVDILLEIYRNNSAGVSDR DSIKA IILDVFAAGTDTTAVVLEWAMTELLRHPEIMKKLQSEVRQVVKDKHNITDDDIEK MHYK AVMKETMRFHPTIPLLVPVARNDVEVMGYDVPVGTMMINAWAIGRDPTSW DEPEKFRP ERFLNSSVDFKGLDFELIPFGAGRRGCPGTTFPMATLEFTLANLMQKFDWELP HECRELD MSERPGVAIRRVIPLLAIGTKM
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	CYP71A8
<b>Protein Names</b>	Recommended name: Cytochrome P450 71A8 EC= 1.14.-.-
<b>Expression Region</b>	1-502
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.