



Recombinant Human T-complex protein 11 homolog (TCP11)

Product Code	CSB-CF855492HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8WWU5
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MPDVKESVPPKYPGDSEGRSCKPETS GPPQEDKSGSEDPPPFLSVTGLTETV NEVSKLSN KIGMNCDYMEKVLPPSSLEGKVKETVHNAFWHDHLKEQLSATPPDFSCALEL LKEIKEI LLSLLLPRQNRRLRIEIEEALDMDLLKQEAHEHGALKVLYLSKYVLNMMALLCAPV RDEAVQ KLENITDPVWLLRGIFQVLGRMKMDMVNYTIQSLQPHLQEHSIQYERAKFQELL NKQPSL LNHTTKWLTQAAGDLTMSPTCPDTS DSSSVAGPSPNEAANNPEPLSPTMVL CQGFLNLL LWDLENEEFPETLLMDRTRLQELKSQLHQLTVMASVLLVASSFSGSVLFGSPQ FVDKLR ITKSLLDFHSRPEEAILTVSEQVSQEIHQSLKNMGLVALSSDNTASLMGQLQNI AKKEN CVCSVIDQRIHLFLKCCLVLGVQRSLLDLPGGLTLIEAELAELGQKFVNLTHHNQ QVFGP YYTEILKTLISPAQALETKVESV
Source	in vitro E.coli expression system
Target Names	TCP11
Protein Names	Recommended name: T-complex protein 11 homolog
Expression Region	1-503
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.