



Recombinant Arabidopsis thaliana WPP domain-interacting protein 2 (WIP2)

Product Code	CSB-CF887751DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9FH18
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MDLESESSVLESVEDNNGLIGDLDKELNSPVETSPLISKGFGLRKWKRRRLRRD LVKDDTS VSMENSKALKRVL SGLVDPNAKQMHLPGPEVRQDSVGSVGSVNSVVGFMVG GESYGNGLA FAAGVDSDNSEDRSSTMHSWDKHRGKVSGGKSVISSGDSSQQRKSSVEKS NKL RGERIK IEKENS HSSMESADSRSSNFVFMQGASYSLS SREQGRRMMDYDDENS DHD AHTSKRKDN VEEEEEETEDYSQGDCVEESQIKSNGSSDNL DPLIVAVNSFQTLQEALQKELQ KFQELGK EEPITSLHDGGESSSCIHAGHEGASEASSSYRFGSEKMGEMELTS LDSEILNLV NNVEHL EIKLEEAKRILEVKETQIRELESTINVSETCNGGTEIGIEDIFQQKVEAEIEYIIFSR SV GNLKRRIKLIEEEKTLGLSKL DKAETKAENLKNQAQDLQNH CVEITEIQEVECLK KRAFK TTRCLLLQLGLLFILYYSLLPEPEIAVPT
Source	in vitro E.coli expression system
Target Names	WIP2
Protein Names	Recommended name: WPP domain-interacting protein 2
Expression Region	1-509
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.