



Recombinant Putative UDP-glucuronosyltransferase ugt-46 (ugt-46)

Product Code	CSB-CF607450CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q10941
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	YKILVFSPATSKSHLISNGRLADELARAGHDVTVLELDFLGISQTTNSVKVAKKR IIDGF QESTNFKNVLHGFSETVMEEPSFTDEIKGWWAYQNVYNDLCAEFLKMDDIFN ELKNAKFD GFFAEQINLCGFGYAHALEIPRHFLISSCPFAAPVYDFTGLPMPTSTVAFAADLS ISPTY TERARNLFVAVLTKLEFTLLNNRLQAHFQHKFGEHFPSLYSVTSDVDVIFVATD EIIDIS TTTLQNIHVGGGLGVDDDDVAEMDNVFASEMSKGKEGVIYFSLGTIANTTKIDSK VMRTVL DIVKKFPDYHFVIRADKYDLSTREYAKSVSNFVSDWLPQPAILHHPRLKLFITH SGYNS IVEAARAGVPLINIPFMFDQNLNSRAVEKKGWGIRRHKKQLLTEPEEIEKAISEII HNKK YSLKAQRIRDLIKSKPLSSSQLLIKTTEWAIKNHGLDEIKFESRGQTTWTYYNLD VIIPV FWLSISLVIPTIFGWYKFSCFGHVEEKKGKSKRD
Source	in vitro E.coli expression system
Target Names	ugt-46
Protein Names	Recommended name: Putative UDP-glucuronosyltransferase ugt-46 Short name= UDPGT 46 EC= 2.4.1.17
Expression Region	18-531
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.