



Recombinant Mouse Fatty acyl-CoA reductase 2 (Far2)

Product Code	CSB-CF745868MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7TNT2
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>MSMIAAFYSNKSILITGATGFLGKVLMEKLFRTSPHLKVIYILVRPKSGQTLQER VFQIL NSKLFKVKVCPNVHEKIRPISADLNQRDFAISKEDVQELLSCTNIIFHCAATV RFDAH LREAVQLNVTATQQLLLMSQMPKLEAFIHISTAFSNCNLSHIDEVIYPCPVEPR KIIDS MEWLDDSIIEEITPKLIGDRPNTYTYTKALGEIVVQQESGNLNVAIVRPSIVGAT WQEPF PGWVDNLNGPSGLIATGKGFRLRSIKATPMAVADVIPVDTVVNLTIAVGWYTAV HRPKST LIYHSTSGNLNPCNWKMGMLQVLATIEKIPFESAFRRPNADFTTSNFTTHYWNT VSHRVP AIIYDFYLRLTGRKPRMLKLMNRLLKTISMLEYFINHSWEWSTNNTEMLLSELS EDQRV FNFDVRQLNWLEYIENYVLGVKKYLLKEDLAGIPKAKQHLRRLRNIHYLFNTALF LIWR LLIARSQMARNVWFFIVSFCYKFISYFRASSTLKV</p>
Source	in vitro E.coli expression system
Target Names	Far2
Protein Names	Recommended name: Fatty acyl-CoA reductase 2 EC= 1.2.1.n2 Alternative name(s): Male sterility domain-containing protein 1
Expression Region	1-515
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.