



Recombinant UPF0392 protein F59C6.8 (F59C6.8)

Product Code	CSB-CF849996CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q93834
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MMLKSCGLIFKGRFVRLFIFIAVCLGFLIAVTILAGLTIFDRQHNHILHDYVARN DDIV VLSTTYEYENSKSFPPNTAVILFNSVQVFHLKYSNLNVVAETMQGNVEVQFKIQP VINTIP FFCKWVPLYLAVGQVPEDHVLLKLSTNKIDGMELSLRTPYETPRKVVACFSPLFL NERWQL LLATVEIYSHYGAFMHFYVRSIITDLFKLIKDNKNTRISPWSAIRIGESRAASPMF DPNT ELEFRNQASAMTDCLLQYKEAAEFIVFPDPDDILVPVLGKNYYEEFTQAFKMFP TAGAVV YNMTQTSIESSMTPALYSPISMLASMKFKGEQKWGKLVVRPERVDSTWIHRSY AIKEGFE QKVMPVDVNAFYHLRIWKFPEVPTFNRSKISNPPFFDPYHLNATKRAIYKISDG LKIQRK FKNRVSQGTMKTIYSRLPKVSLYYPLIEVCYNRIFYSMKDIGTCRGPEYCNIPAF PGLRC TNVASEFVYKSYRNIYIHQLISTDFEEGDNGCTL
Source	in vitro E.coli expression system
Target Names	F59C6.8
Protein Names	Recommended name: UPF0392 protein F59C6.8
Expression Region	1-515
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.