



Recombinant Escherichia coli Uncharacterized protein YlaB (ylaB)

Product Code	CSB-CF302054ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77473
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	<p>MRTRHLVGLISGVLILSVLLPVGLSIWLAHQVETSFIEELDTYSSRVAIRANKV ATQ GK DALQELERWQGAACSEAHLMEMRRVSYSYRYIQEVAYIDNNVPQCSSLEHES PPDTFPEP GKISKDGYRVWLTSNDLGIIRYVMAMGTAHYVVMIDPASFIDVIPYSSWQIDA AIIGNA HNVVITSSDEIAQGIITRLQKTPGEHIENNGIYDILPLPEMNISIITWASTKMLQKG WH RQVFIWLPLGLVIGLLAAMFVLRILRRIQSPHHRLQDAIENRDICVHYQPIVSLAN GKIV GAEALARWPQTDGSWLSPDSFIPLAQQTGLSEPLLLIIRSVFEDMGDWLRQH PQQHISI NLESPVLTSEKIPQLLRDMINHYQVNPRQIALELTEREFADPKTSAPIISRYREA GHEIY LDDFGTGYSSLSYLQDLVDILKIDKSFVDALEYKNVTPHIIEMAKTLKLMVAE GIETS KQEEWLRQHGCVHYGQGWLYSKALPKEDFLRWAEQHL</p>
Source	in vitro E.coli expression system
Target Names	pdeB
Protein Names	Recommended name: Uncharacterized protein YlaB
Expression Region	1-516
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.