



# Recombinant Putative UDP-glucuronosyltransferase ugt-58 (ugt-58)

<b>Product Code</b>	CSB-CF630288CXY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q20086
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Sequence</b>	<p>GDVLFIPSTLYPVHAQTMSVLARELAERGHKVTWLEIGPERSDIVLHQQVSREF  WPAQFG  DRTLQEIQFKNHTSHSELWNPQYTNENEQTTGWLASIRLCDSVLSRSKSKFD  KMVEKQF  STVIVDDLYNSCGVLMAGLKKSVMYIYWSMTGLRTEAWANQSPSPPSYLPVA  GTGLTDDL  TFSQRLYNVASYFKQLYIHQHIVQPRIDAVFQKYYPGVESTFEIERNASINFVNT  PPIFD  FSRPYMPRVNFIGAIQCRKPKELPKEFTSWISAYPDGFVVLSTGFTVQWNKSP  EHIRQAY  LDTFKYLPNLLFIWQTGLPNSSNLPSNLLVKPWLPLQDLLGHQKCRCHVSHGG  LNSVIES  VYHGVPVVGVPPLTSRGYDNLRLRITARDSGVMVEKSEFSGEILTAINEVIENEK  YKKEML  IFQDMVIDVPYTELYHAAFVVEFIERHQEVPHARSGADHLNFLQYFLVDVIAFF  FFVIFC TLSIIYYTIRTVFRTIRSVVNGVIRGVPKIVSRGKNN</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	ugt-58
<b>Protein Names</b>	Recommended name: Putative UDP-glucuronosyltransferase ugt-58 Short name= UDPGT 58 EC= 2.4.1.17
<b>Expression Region</b>	17-533
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.