



Recombinant Human Motile sperm domain-containing protein 2 (MOSPD2)

Product Code	CSB-CF839916HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NHP6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MAENHAQNKAKLISETRRRFEAEYVTDKSDKYDARDVERLQQDDNWWVESYLS WRHNIVDE TLKMLDESFQWRKEISVNDLNESSIPRWLLEIGVIYLHGYDKEGNKLFWIRVKY HVKDQK TILDKKLIAFWLERYAKRENGKPVTVMFDLSETGINSIDMDFVRFIINCCKVYYP KYLS KIVIFDMPWLMNAAFVKTWLGPEAVSLLKFTSKNEVQDYVSVEYLPPHMG TDPFKYS YPPLVDDDFQTPCENGPISEDETSSKEDIESDGKETLETISNEEQTPLLKIN PTEST SKAEENEKVDKVKAFKKPLSVFKGPLLHISPAEELYFGSTESGEKKTIVLTVN TKNIV AFKVRTTAPEKYRVKPSNSSCDPGASVDIVVSPHGGLTVSAQDRFLIMAAEME QSSGTGP AELTQFWKEVPRNKVMEHRLRCHTVESKPNLTLKDNAFNMSDKTSEDICL QLSRLLES NRKLEDQVQRCIWQQLLLSLTMLLLAFVTSFFYLLYS</p>
Source	in vitro E.coli expression system
Target Names	MOSPD2
Protein Names	Recommended name: Motile sperm domain-containing protein 2
Expression Region	1-518
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.