



Recombinant UPF0392 protein C35A5.5 (C35A5.5)

Product Code	CSB-CF613830CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q18473
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MTKIKRQINNFDFRRIKILLKKWNVVIYFVLILICFYFIPIYFPNNDKMNLWLSSIK YY LTYPIYNESLTKTDAYIINTYYPKSSSLGENAIGMILLMNRNTQRDMTKYRMKLI ASNS SHQSVIVTPKFLEESYSSCPYINMVAMVNTLPNLNKLEIFDGERKMEIPFQMKG TTAPAS VIICIS PQFVAEQWQLFVAHAHVARKFGGHLHMYVTSIIDTFFDLVQEYERLGY VTIDYW MRLKLANSSVDSVEPNLHSELRNQAGA QSDCLYQYKEAAAFITFFDLDDIFIPR GYDSYF DEFSALYELHPNILTFQYTKRETMVYNKAKIEDINFEELFGHTWVNEEDYGKV MTKPGN LNSMWIHESWNIPTNRHHVSKSNYIIHMQKPVDPDGTDPVSYRMSNFEMLES MQLNASTL IPIQDDLRLV LKSSNISMTAEDLPNKTYYPFIYRCYYEKFYKKPKKDSCPNGEG CLIPQ RTGTNCIHSADDFKSGPEMWPITYHYHVNSKWSRTKGCHA
Source	in vitro E.coli expression system
Target Names	C35A5.5
Protein Names	Recommended name: UPF0392 protein C35A5.5
Expression Region	1-520
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.