



Recombinant *Saccharomyces cerevisiae* Formation of crista junctions protein 1 (FCJ1)

Product Code	CSB-CF451447SVP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3LRA0
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Sequence	<p>ASINTGTTVASKKASHKFRNTFWTIALSATAFYAGGIIYSQKNDKFGDFFSNNV PFAEDL LETYEHYHDRPTLFLEDSWDGLKAKSNDLLSGLTGSSQTRRSNRENIEVKKILS LEPLNI ETENSDPQLKEIIGSLNDLINSNLSNLSIPSEFNSIKKSNQNMLTNLSQLNETL KEAL SNYMIQRTSEVITELNTQYENSKREFEKNLQKNLLQEVDEFKENLTKQKDKLE EKLKAN EELLQAKHANEVGLLSITQVKEFNKIIKDKIEKERNGLAHLEEINSEVNDLSKSI DRSS KILSKNEALVQLTFQVDEIKSRINNNLPDVNIDKELSRLKLLSNLLSTFNKKSCC DDGD CCSCKKGNKNEGKEGKISCKCKPKTNPPSLLSVALDELESTCSGKKILSNEQIY NRWNLL ADDFKTASLLPPNSGILGQLTAKVFSFLFTKTGNPSNATDFDSVYARVGDNLR VSNLND AVEEVVSLKGWPHKVCESWIEDARRKLEVQRLVEILDCEIRTL</p>
Source	in vitro E.coli expression system
Target Names	MIC60
Protein Names	Recommended name: Formation of crista junctions protein 1 Alternative name(s): Altered inheritance of mitochondria protein 28 Found in mitochondrial proteome protein 13 Mitofilin
Expression Region	17-539
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.