



# Recombinant Bovine Chloride channel CLIC-like protein 1 (CLCC1)

<b>Product Code</b>	CSB-CF638074BO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q1LZF8
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Sequence</b>	<p>HDDEWIDPTDMLNYDAASGRMRKSQVKYGISEKEEVNPDLSANELSECYNR  LDSLTYKI  DECEKQKRKDYESQSNPVFRRYLNKILTIETKKLGLPDENKHMHYDAEIIILKRQ  TLLEIQ  KFLSGEDWKPGALDDALSDILINFKFHDFETWKWRFEEFFGVDPYNVFMVLLC  LLCIVAL  VATELWTYVRWYTQLKRVFFISFLISLGNWWMYLYKLAFQAQAEVAKMEPLN  NVCAEKM  NWSGSLWEWLRSSWTYKDDPCQKYELLVNPIWLVPTKALAVTFTNFVTE  PLKHVGKG  AGEFIKALMKEIPVLLHIPVLIIMALAVLSFCYGAGKSVNMLRHVGGPEREAPQA  LQAGE  RRRQQKIDYRPHGGAGDADFYYRGQISPIEQGPNDNTYEGRRDVLRRERDVGL  RFQTGNKS  PEVLRPFDLQEAAREHPKVVPGLKSPNLESKPREMGEIPGESTPTESSTESS  QPAKPVS GQKVSEGVEGCPAVEKAQLRRTDAAGGPEEGSTCSPASTAVEVCG</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	CLCC1
<b>Protein Names</b>	Recommended name: Chloride channel CLIC-like protein 1
<b>Expression Region</b>	19-542
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.