



Recombinant Danio rerio DDB1- and CUL4-associated factor 17 (dcaf17)

Product Code	CSB-CF003842DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8WHR0
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>MCALPSPRAQRPPACTRRTSKNSCALLASRYNGSCSNDFGRLLASNLKILRNIL KDDTEF VKVWSKTSKSLISYESGRIYFDNYRCCYSSLLPEPELLYELPRTPKTEKIEDALL CQCPL ENVLPNASEQKSCLLTLTANNWLYLLSADTGETLQRVYLSSRFKFSRSLGWDS SQETFYV KSVQCKQTALERQAGVDNINILMRVAAFQVFPLQLVGMLEINKKVFGKTAVDVL LSQGVLA VSHSSKLVRLYSFEYIVNKFRTTELILGQQCELNNARGIVGDAPYGVPVNICIHE CPPVL FEMTYFENGVQIGGHPWHYIYTPNHKRHRGTHHVCSITDGALAKNGVQDMKC DSLEVDWI FFHPDDSGRIIHAGPSTINILKIMAETGCDWKYEIITDFSITAARDSNASQVIVTSS GRT VKRRFQMLDDDPAQETFRMVKYEDELDLLAVVDITHAEDEGQARLRLHDNKT GALMKRVP LEEPWDVTYSHEVYFDRDTHIHTVQEKNNSFCCHVYKMKRPASDQP</p>
Source	in vitro E.coli expression system
Target Names	dcaf17
Protein Names	Recommended name: DDB1- and CUL4-associated factor 17
Expression Region	1-526
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.