



Recombinant Escherichia coli Uncharacterized membrane protein YjcC (yjcC)

Product Code	CSB-CF330190ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32701
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MSHRARHQLLALPGIIFLVLPFIILSLWIAFLWAKSEVNNQLRTFAQLALDKSELV IRQA DLVSDAAERYQQGVCTPAHQKRMLNIIRGYLYINELIYARDNHFLCSSLIAPVN GYTIAP ADYKREPNSIYYYYRDTPPFFSGYKMTYMQRGNVAVINPLFWSEVMSDDPTL QWGVYDTV TKTFFSLSKEASAATFSPLIHLKDLTVQRNGYLYATVYSTKRPIAAIVATSYQRLLI THFY NHLIFALPAGILGSLVLLLLWLRIRQNYLSPKRKLQRALEKHQLCLYYQPIIDIKTE KCI GAEALLRWPGEQQQIMNPAEFIPLAEKEGMIEQITDYVIDNVFRDLGDYLATHA DRYVSI NLSASDFHTSRLIARINQKTEQYAVRPQQIKFEVTEHAFLDVKMTPILAFRQA GYEVA IDDFGIGYSNLHNLKSLNVDILKIDKSFVETLTTHKTSHLIAEHIIELAHSLGLKTIA EG VETEEQVNWLRKRGVRYCQGWFFAKAMPPQVFMQWMEQLPARELTRGQ
Source	in vitro E.coli expression system
Target Names	pdeC
Protein Names	Recommended name: Uncharacterized membrane protein YjcC
Expression Region	1-528
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.