



Recombinant *Saccharomyces cerevisiae* Autophagy-related protein 32 (ATG32)

Product Code	CSB-CF408155STA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6ZVD0
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Sequence	<p>MVLEYQQREGKGGSSSKSMPPDSSSTTIHTCSEAQTGEDKGLLDPHLSVLELLS KTGHSPS PMGQNLVTSIDISGNHNVNDSISGSWQAIQPLDLGASFIPERCSSQTTNGSILS SSDTSE EEQELLQAPAAADIINI KQGQEGANVVSPSHPFKQLQKIISLPLPGKEKTPFNEQ DDDGD EDEAFEEDSVTITKSLTSSSTNSFVMPKLSLTQKNPVFRLILGRTGSSFYQSIPK EYQSL FELPKYHDSATFPQYTGIVIFQELREMVSLNLRIVQYSQGKPVIPICQPGQVIQV KNLV KSFLRNKLVKLLFPPVVVTNKRDLKMFQRLQDLSLEYGEDVNEEDNDDEAIH TKRSYC RNKKAENSKKKSPKSNKKPKRKKQKFFTSWFTWGISITIGISFGCCVTYFVTAA YEHQTV KSLSLRPSILASLLSLDSSSDTINTPATASPSSTEQFLWFDKGTQINFHSDGFI MKSLT IIKETWGMNTFVLHALSKPLKFLLENLKSSEFSIDESNRILALGYILL</p>
Source	in vitro E.coli expression system
Target Names	ATG32
Protein Names	Recommended name: Autophagy-related protein 32 Alternative name(s): Extracellular mutant protein 37
Expression Region	1-529
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.