



Recombinant Protein ST7 homolog (F11A10.5)

Product Code	CSB-CF618296CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q19337
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MACSWTFLWLLWIALVAVLLFALRGPLKISESLESVTATSYFNNLTPKFYVALT GTSSLV SGIILIFEWWYFKNNAGIDAGDEEGSDNDESIENTKTVPECKVWRNPMALFRA AEYNRFR KETNSEPLTYDDMNLQAQDHQSLFMCDEDQGRAEYEIMQVAWRERERESEERIQ TARAALAI NPECASALVLLAEEEESETVSAENLLRRALRAIESTLNSYSNNQIASYAQNGDA VRKRDL TIQTYIKRRLAMCARKQGRLREAIGFRDLSRDQSLSTLLSVQDNLEACLEVQ AYADVQ NLLVRYDGYGAPCSYELREPRSAAMSYSALLKVRVAENFRCAADSSIRRL SSAEQTA IEALTRAMEFNPHVPPYLLELRSMIMPPEHFLKRGDSEALAYAFFHIQHWKRID GALQLL SIVWKDFVPKVSKDKNAFSSQLESADKELLPSWHEQSVFPKTEGLMMLLQTF ICLAICI LAVLAQQFPASSGEIFRTAATIGMQFYENSVYTVSQWAPGNIIPYLASKQVPVP EL
Source	in vitro E.coli expression system
Target Names	F11A10.5
Protein Names	Recommended name: Protein ST7 homolog
Expression Region	1-536
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.