



Recombinant Human Atlastin-3 (ATL3)

Product Code	CSB-CF728503HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6DD88
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLSPQRVAAAASRGADDAMESSKPGPVQVVLVQKDQHSFELDEKALASILLQ DHIRDLDV VVVSVAGAFRKGKSFILDFMLRYLYSQKESGHSNWLGDPEEPLTGFSWRGGG DPETTGIQ IWSEVFTVEKPGGKKA VVLMMDTQGAFDSQSTVKDCATIFALSTMTSSVQIYNL SQNIQE DDLQQLQLFTEYGRLAMDEIFQKPFQTLMFLVRDWSFPYEYSYGLQGGMAFL DKRLQVKE HQHEEIQNVRNHIHSCFSDVTCFLLPHPGLQVATSPDFDGLKLDIAGEFKEQL QALIPYV LNPSKLMEKEINGSKVTTCRGLLEYFKAYIKIYQGEDLPHPKSMLQATAEANNLA AAASAK DIYYNMEEVCGGEKPYLSPDILEEKHCEFKQLALDHFKTKKMGKDFSFY QQELEEE IKELYENFCKHNGSKNVFSTFRTPAVLFTGIVALYIASGLTGFIVLEVAQLFNC MVGLL LIALLTWGYIRYSGQYRELGGAIIDFGAAYVLEQASSHIGNSTQATVRDAVVGRP SMDKKA Q
Source	in vitro E.coli expression system
Target Names	ATL3
Protein Names	Recommended name: Atlastin-3 EC= 3.6.5.-
Expression Region	1-541
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.