



Recombinant Human Protein FAM134A (FAM134A)

Product Code	CSB-CF850856HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NC44
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MASGGGGGNTGAGGGPGMGLSLGLGLSLGLMSEATSEAEAAAATAEAVGR LATTLWLRL RGWEAVLAAAQRLLVWEKPLHSLVTAALNGLFWLLSSSSLRPFFLLSVSLLA YFLDLW QPRFLPDVSASSPEEPHSDSEGAGSGARPHLLSVPELCRYLAESWLTQIHLQ ELLQYKR QNPAQFCVRVCSGCAVLAVLGHYVPGIMISYIVLLSILLWPLVVYHELIQRMYTR LEPLL MQLDYSMKAEANALHHKHDKRKRQGKNAPPGGDEPLAETESEAEELAGFS PVVDVKKTA LALAITDSELSDEEASILESGGFSVSRATTPQLTDVSEDLDQQSLPSEPEETLS RDLGEG EEGELAPPEDLLGRPQALSRQALDSEEEEEEDVAAKETLLRLSSPLHFVNTHFN GAGSPPD GVKCSPPGPVETLSPETVSGGLTALPGTLPPLCLVGS DPAPSPSILPPVPQD SPQPLPA PEEEEEALTTEDFELLDQGELEQLNAELGLEPETPPKPPDAPPLGPDHSLVQSD QEAQAV AEP</p>
Source	in vitro E.coli expression system
Target Names	RETREG2
Protein Names	Recommended name: Protein FAM134A
Expression Region	1-543
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.