



Recombinant Human Reticulon-2 (RTN2)

Product Code	CSB-CF020569HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75298
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>MGQVLPVFAHCKEAPSTASSTPDSTEGGNDDSDFRELHTAREFSEEDDEEETT SQDWGTPR ELTFSYIAFDGVVGSRRDSTARRPRPQGRSVSEPRDQHPQPSLGDSLESI PSLSQSPE PGRRGDPDTAPPSEPLEDLRLRLDHLGWVARGTGSGEDSSTSSSTPLEDEE PQEPNRLE TGEAGEELDLRLRLAQSSPEVLTPQLSPGSGTPQAGTPSPSRSRDSNSGPE EPLLEEEE KQWGPLEREPVRGQCLDSTDQLEFTVEPRLLGTAMEWLKTSLLLAVYKTVPIIL ELSPLW TAIGWVQRGPTPTPVLRVLLKWAKSPRSSGVPSLSLGADMGSKVADLLYWK DTRTSGVV FTGLMVSLCLLHFSIVSVAHLALLLLCGTISLRVYRKVLQAVHRGDGANPFQ AYLDVD LTLTREQTERLSHQITSRVVSAATQLRHFFLVEDLVDSLKLALLFYILTFVGAIFN GLTL LILGVIPLFTIPLLYRQHQAIQIDQYVGLVTNQLSHIKAKIRAKIPGTGALASAAAA VSGS KAKAE</p>
Source	in vitro E.coli expression system
Target Names	RTN2
Protein Names	Recommended name: Reticulon-2 Alternative name(s): Neuroendocrine-specific protein-like 1 Short name= NSP-like protein 1 Neuroendocrine-specific protein-like I Short name= NSP-like protein I Short name= NSPLI
Expression Region	1-545
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.