



Recombinant Escherichia coli O127:H6 Probable ubiquinone biosynthesis protein UbiB (ubiB)

Product Code	CSB-CF484801EOB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7UNG5
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Sequence	<p>MTPGEVRRLYFIIRTFSLYGLDELIPKMRITLPLRLWRYSLFWMPNRHKDKPLG ERLRLA LQELGPVWIKFGQMLSTRDLFPPHIADQLALLQDKVAPFDGKLAQQIEAAM GGLPVEA WFDDFEIKPLASASIAQVHTARLKSNGKEVVIVIRPDILPVIKADLKLILYRLARW VPRL LPDGRRLRPTEVVREYEKTLIDELNLLRESANAIQLRRNFEDSPMLYIPEVYPD YCSEGM MVMERIYGIPVSDVAALEKNGTNMKLLAERGVQVFFFTQVFRDSFFHADMHPG NIFVSYEH PENPKYIGIDCGIVGSLNKEDKRYLAENFIAFFNRDYRKVAELHVDSGWVPPDT NVEEFE FAIRTVCEPIFEKPLAEISFGHVLLNLFNTARRFNMEVQPQLVLLQKTLLEYEGV GRQLY PQLDLWKTAKPFLESWIKDQVGIPALVRAFKEKAPFWVEKMPPELPELVYDSL RQGKYLQH SVDKIARELQSNHVRQGQSRVFLGIGATLVLSGTFLVSRPEWGLMPGWLMA GGLIAWFV GWRKTR</p>
Source	in vitro E.coli expression system
Target Names	ubiB
Protein Names	Recommended name: Probable ubiquinone biosynthesis protein UbiB
Expression Region	1-546
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.