



Recombinant Salmonella paratyphi C Probable ubiquinone biosynthesis protein UbiB (ubiB)

Product Code	CSB-CF507171SWU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C0Q3E3
Product Type	Transmembrane Protein
Immunogen Species	Salmonella paratyphi C (strain RKS4594)
Sequence	<p>MTPGEVRRLYFIIRTFSLYGLDELIPRMRLTLPLRLWRYSLFWMPNRHKDKLLG ERLRLA LQELGPVWIKFGQMLSTRDLFPPQIADQLALLQDKVAPFDGRLAKAQIEEAM GGLPVEA WFDDFDIQPLASASIAQVHTARLKSNGKEVVIVIRPDILPVIQADLKLIIYRLARW VPRL LPDGRRLRPTEVVREYEKTLIDELNLLRESANAIQLRRNFENSPMLYIPEVYSD YCSQNM MVMERIYGIPVSDVAALEKNGTNMKLLAERGVKVFFTQVFRDSFFHADMPHG NIFVSHEH PENPQYIGIDCGIVGSLNKEDKRYLAENFIAFFNRDYRKVAELHVDSGWVPPDT NVEDFE FAIRTVCEPIFEKPLAEISFGHVLLNLFNTARRFNMEVQPQLVLLQKTLLEYEGV GRQLY PQLDLWKTAKPFLESWIKDQVGIPALTRALKEKAPFWVEKMPEIPELVYDSLRLQ GKYLQH SVDKIARELQVNHVRQSQRSRYLLGIGATLLLSGSFLLVNRPEWGLMPGWLMV GGVVVWLV GWRKTR</p>
Source	in vitro E.coli expression system
Target Names	ubiB
Protein Names	Recommended name: Probable ubiquinone biosynthesis protein UbiB
Expression Region	1-546
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.