



Recombinant Escherichia coli Paraquat-inducible protein B (pqjB)

Product Code	CSB-CF331911ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P43671
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	<p>MESNNGEAKIQKVKNWSPVWIFPIVTALIGAWVLFYHYSHQGPEVTLITANAEG IEGGKT TIKRSVDVGVVESATLADDLTHVEIKARLNSGMEKLLHKDTVFWVVKPQIGRE GISGLG TLLSGVYIELQPGAKGSKMDKYDLLDSPPLAPPAKIRVILDSKKAGQLSPGD PVLFRG YRVGSVETSTFDTQKRNIYQLFINAPYDRLVTNNVRFWKDSGIAVDLTSAGM RVEMGSL TTLLSGGVSFVDVPEGLDLGQPVAPKTAFLYDDQKSIQDSLYTDHIDYLMFFKD SVRGLQ PGAPVEFRGIRLGTVSKVPPFAPNMRQTFNDDYRIPVLIRIEPERLKMQLGENA DVVEHL GELLKRGLRGSLKTGNLVTGALYVDLDFYPNTPAITGIREFNGYQIIPTVSGGLA QIQQR LMEALDKINKLPLNPMIEQATSTLSESQRMTMKNLQTTLD SMNKILASQSMQQLP TDMQST LRELNRSMQGFQPGSAAYNKMVADMQRDLQVLRQLPVLKTLNEKSNALVFE AKDKKDPE PKRAKQ</p>
Source	in vitro E.coli expression system
Target Names	pqjB
Protein Names	Recommended name: Paraquat-inducible protein B
Expression Region	1-546
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.