



Recombinant UPF0392 protein R07B7.12 (R07B7.12)

Product Code	CSB-CF630354CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q21802
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MGNFSFTMFQILNLKMNLSPYCNRSTGRICSKIRLYHRSFRSNPSLQMCFIIVF LLFIF SLIMFMKLNQNSYSYGKEILGYLYETAIDQLTYNHHTAFITSAYYYRNSKSLGENA VAMIV AMSQRTEFKHLENHEITIVGSRGYQNLKTKASITVETAEQMSCEYVMALIQQGITLE IPQKL EIESAGTRVQIPFREPRKNSHSPVIICISPPQFVAEKWQLFLMNIHVIRRYGGHMH IYITS MVEKLFNVLKIYEDMEALTIDYWIRMKLKKTSSPVADIMKNVEWRHQAGAQT DCLLYKE VAEFIAFFDIDDILVPNFHSHYHQEFSSHFNAFPSYHSIFYGKRDFVFEKISSIED FSFR HLFSNMKIQEETGYGKSIVNPLKYNSTWIIHSMKLPKMLKIMNTEIIHIKNILD SELN KDAPIHLPPIYGTETESVIREMDLKTLDLDFQTVYQNPIYREAALKMIDYNFYTP I VFNC YNESFYHPYFVEKKDFSQICPNADNCQLPQREDIKCIHSDGEYVSGPEMYPITF HYAVHP FWSNDIGCYQ
Source	in vitro E.coli expression system
Target Names	R07B7.12
Protein Names	Recommended name: UPF0392 protein R07B7.12
Expression Region	1-550
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.