



Recombinant Pongo abelii Abhydrolase domain-containing protein 16A (ABHD16A)

Product Code	CSB-CF712854PYX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5R6S0
Product Type	Transmembrane Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sequence	AKLLSCVLGPRLYKIYRERDSEAPASVPETPTAVTAPHSSSWDTYYQPRALE KHADSIL ALASVFWISISYSSPFAFFYLRYKGYLSLSKVVVPSHYAGTLLLLLAGVACLRGI GRWTN PQYRQFITILEATHRNQSSSENKRQLANYNFDFRSWPVDFHWEEPSSSRKESRG GPSRQGVA LLRPEPLHRGTADTLLNRVKKLPCQITSYLVAHTLGRMLYPGSVYLLQKALMP VLLQGQ ARLVEECNGRRAKLLACDGNEIDTMFVDRRGTAEPQGQKLVICCEGNAGFYE VGCVSTPL EAGYSVLGWNHPGFAGSTGVFPQNEANAMDVVVQFAIHR LGFQPQDIIIIYAW SIGGFTA TWAAMSYPDVS AVILDASFDDL VPLALKVMPDSWRGLVTRTRVRQH LNLNNAE QLCRYQGP VLLIRRTKDEIITTTVPEDIMSNRGNDLLLKLLQHRYPRVMAEEGLRVVRQWLE ASSQLE EASIYSRWEVEEDWCLSVLRSYQAEHGPDFPWSVGEDMSADGRRQLALFLA RKHLHNFEA THCTPLPAQNFQMPWHL
Source	in vitro E.coli expression system
Target Names	ABHD16A
Protein Names	Recommended name: Abhydrolase domain-containing protein 16A EC= 3.-.- Alternative name(s): HLA-B-associated transcript 5
Expression Region	2-558
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.