



Recombinant Human Protein FAM55B (FAM55B)

Product Code	CSB-CF853402HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96DL1
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MVEKILIHRLTLFPNAIARKLLLMLTFILIFWIIYLASKDHTKFSFNLENHII LNQGN I FKKYSHSETPLCPAVSPKETELRIKDIMEKLDQQIPPRPFTHVNTTTSATHSTAT ILNPQ DTYCRGDQLDILLEVRDHLGHRKQYGGDFLRARMYSTALMAGASGKVTD FNN GTYLVSFT LFWEGQVSLSLLLIHPSEGV SALWRARNQGCDRIIFTGLFANRSSNVFTECGLT LNTNAE LCQYMDDRDQEA FYCVRPQHMPCEALTHMTTRTRNISYLSKEEWRLFHRSNI GVEMMKNF TPIEVIPCNKSENIKKNCQIGMKTPFPSGYTLKKMWITAFCKQIKFNETKNINDC LERKL IYLMGDSTLHQWIYYLQKAVKTLKYFDHHGAGIFKTHVLLDVERHILIQWKKHG HPFVTK KLFSVKDENYIPREIDQVAGDKNTAIVITLGQHFRFPINIFIRRAINIQA IERLFL RS PETKVILKTENTREIEQNAEMFSDFHGYIQNLIIRDIFVDLNVGIIDAWDMTIAYCT NNA HPPDYVIQNQIGMFLNYIC
Source	in vitro E.coli expression system
Target Names	NXPE2
Protein Names	Recommended name: Protein FAM55B
Expression Region	1-559
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.