



Recombinant Salmonella typhimurium Flagellar M-ring protein (fliF)

Product Code	CSB-CF323397SXB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15928
Product Type	Transmembrane Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Sequence	<p>SATASTATQPKPLEWLNRLRANPRIPLIVAGSAAVAIVVAMVLWAKTPDYRTLFSNLSDDGGAIVAQLTQMNIPYRFANGSGAIEVPADKVHELRLRLAQQGLPKGGAVGFELLDQEKFGISQFSEQVNYQRALEGEARTIETLGPVKSARVHLAMPKPSLFVREQKSPSASVTVTLEPGRALDEGQISAVVHLVSSAVAGLPPGNVTLVDQSGHLLTQSNTSGRDLNDAQLKFANDVESRIQRRIEAILSPIVGNGNVHAQVTAQLDFANKEQTEEHYSPNGDASKATLRSRQLNIS</p> <p>EQVGAGYPPGGVPGALSNDPAPPNEAPIATPPTNQQNAQNTPTSTSTNSNSAGPRSTQRNETSNYEVDRTIRHTKMNVDIERLSVAVVNYKTLADGKPLPLTADQMKQIEDLTREAMGFSDKRGDTLNVVNSPFSAVDNTGGELPFWQQQSFIDQLLAAGRWLLVVLVAVWILWRKAVRPQLTRRVEEAKAAQEQAVRQETEEAVEVRLSKDEQLQRRANQRLGAEVMSQRIREMSDNDPRVVALVIRQWMSNDHE</p>
Source	in vitro E.coli expression system
Target Names	fliF
Protein Names	Recommended name: Flagellar M-ring protein
Expression Region	2-560
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.