



Recombinant Human Transmembrane protein 151B (TMEM151B)

Product Code	CSB-CF023727HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8IW70
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSPPGSAAGESAAGGGGGGGPGVSEELTAAAAAAADEGPAREEQRPPIQP SFTKSLCRE SHWKCLLLSLLMYGCLGAVAWCHVTTVTRLTFSSAYQGNSLMYHDSPCSNGY VYIPLAFL LMLYAVYLVECWHCQARHELQHRVDVSSVRERVGRMQQATPCIIWWKAISYH YVRRTQVT RYRNGDAYTTTQVYHERVNTHTVAEAEFDYARCGVRDVSKTLVGLEGAPATRL RFTKCFSE ASVEAENAYLCQRRFFAENEGLDDYMEAREGMHLKNVDFREFMVAFPDPA RPPWYACSS AFWAAALLTLWPLRVLAEYRTAYAHYHVEKLFGLGPGSASSAGGGLSPSD ELLPLTH RLPRVNTVDSTELWHIRSNQQLVPSYSEAVLMDLAGLGTRCGGAGGGYAPS CRYGGVGG PGAAGVAPYRRSCEHCQRAVSSSSIFSRSAICASPRAGPGPGGGAGCGGS RFSLGRLY GSRRSCLWRSRSGSVNEASCPTQTRLSSQASMGGDEDDDEEEAGPPPPYH DALYFPVLI VHRQEGCLGHSHRPLHRHGSCVETSL
Source	in vitro E.coli expression system
Target Names	TMEM151B
Protein Names	Recommended name: Transmembrane protein 151B Alternative name(s): Transmembrane protein 193
Expression Region	1-566
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.