



# Recombinant *Saccharomyces cerevisiae* Uncharacterized protein SLP1 (SLP1)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-CF619573SVG  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | Q12232   |
| <b>Product Type</b>      | Transmembrane Protein  |
| <b>Immunogen Species</b> | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)   |
| <b>Sequence</b>          | LDRNKTAQNAKIGLHDTTVITTGSTTNVQKEHSSPLSTGSLRTHDFRQASKVDI<br>RQADIR<br>ENGERKEQDALTQPATPRNPGDSSNSFLSFDEWKKVKSKEHSSGPERHLSR<br>VREVPDPSC<br>YKEKECIGEELEIDLGFLTnkNEWSEREENQKGFNEEKDIEKVYKKKFNYASLD<br>CAATIV<br>KSNPEAIGATSTLIESKDKYLLNPCSAPQQFIVIELCEDILVEEIEIANYEFFSSTF<br>KRF<br>RVSVSDRIPMVKNEWTLGFEFARNRSRELQKFQIHNPQIWASYLKIEILSHYEDE<br>FYCPI<br>SLIKVYGKSMDEFKIDQLKAQEDKEQSIGTNNINNLNEQNIQDRCNNIETRLE<br>TPNTSN<br>LSDLAGALSCTSKLIPLKFDEFFKVLNASFCPSKQMISSSSSSAVPVIPEESIFKN<br>IMKR<br>LSQLETNSSLTVSYIEEQSKLLSKSFEQLEMAHEAKFSLVTFIFNETMMSNLDLL<br>NNFAN<br>QLKDQSLRILEEQKLENDKFTNRHLLHLERLEKEVSFQRRIVYASFFAFVGLISY<br>LLITR ELYFEDFEESKNGAIEKADIVQQAIR |
| <b>Source</b>            | in vitro E.coli expression system  |
| <b>Target Names</b>      | SLP1   |
| <b>Protein Names</b>     | Recommended name: Uncharacterized protein SLP1 Alternative name(s): SUN-like protein 1   |
| <b>Expression Region</b> | 22-587   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full Length of Mature Protein  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |