



Recombinant Rat WSC domain-containing protein 1 (Wscd1)

Product Code	CSB-CF700534RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q505J3
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>MAKPPFRLQKFLRRTQFLLLFLTAAYLMTGSLLLLQRARVALPQALRAPGSLQA LPVATV ALGVGLLDSRSLHDPHSSPDLLLDVDTLRSPLARVPPGIRWPRRNRSSLRRR WLHHLTSD PQGPPTLGPEASGPASHNQGNYLGCFFSEEGQERTLKGAVFYDLRKMTVSHC QEACAERSY VYAGLEAGAECYCGNRLPATRVSLKECNQECKGEKGSMEGALHRLSVYSVG LQQSGAKKR WTATYRGCFPLPENITRTFSSSMTQANVTVETCSGFCSQKEFPLAILRGWACY CAYTPPQ FSLRDAVDGALCSQGTEAQGLPGYCEVYQTPVQDTRCTDRKFLPNKSKVFVA LSSFPGAG NTWARHLIEHATGFYTGSIYFDGTLYNKGFKEKDHWSRRTICVKTHESGR REIEMFDS AILLIRNPYRSLVAEFNRKCAGHLGYAPDRNWKSKEWPDFVNSYASWWSSHV LDWLKYGK RLLVVHYEELRHSLVPTLREMVAFLNVSVEERLLCVENKESFRRRGRPHPH DQEPFTP EMKDLINGYIRTVDQALRDHNWAGLPREYVPR</p>
Source	in vitro E.coli expression system
Target Names	Wscd1
Protein Names	Recommended name: WSC domain-containing protein 1
Expression Region	1-572
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.