



Recombinant Danio rerio WSC domain-containing protein 2 (wscd2)

Product Code	CSB-CF026156DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2BGL3
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>MARSLLLKIHRYFRRKPVRRFFSFILLYLTAGSLVFLHSGFSSDSSTAGIAGSGSDP LLSEG RGGTGGGSATSEGLGLLGRVFKETRRTPRRFGPPWMKESRGQDAPEWAGR GFDHTSSWSH GAKGRRTKEMDDGRAKYIGCYVDDTQKRALRGVSFLDYKKMTVFRQCNDCA ERGYLYAGL EFGAECYCGHKIQAPNVSESECNMECKGEKSNLCCGGPNRLSIYRLELSQESA RRYGSAIF KGCRRPDNVTALPASAVMQNMSVDKCVDMCTEKEFSLAALAGDKCHCGF PTPLFNLHE HEDEELCLHRCTGEDFESCGNRDFFVYQTQVQDNRCMDRRFLPTRSKHLM ALASFPGAG NTWARHLIELATGYTGSYYFDGSLYNKGFKGERDHWRSGR TICIKTHESGKK EIETFDA SILMIRNPYKALMAEFNRKYGGHIGFASQAHWRGKEWPEFVKNYAPWWASHT LDWLKYGK KVQVVHFEDLKRDLFSQLKGMVIFLGLVSEDRLLCVEGQKDGNFKRSGLRKL EYDPYTV EMRATIDRLIKTVDMELRKRNLGVPDEYRPR</p>
Source	in vitro E.coli expression system
Target Names	wscd2
Protein Names	Recommended name: WSC domain-containing protein 2
Expression Region	1-572
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.