



Recombinant Oncorhynchus mykiss Complement component C9 (c9)

Product Code	CSB-CF361918OEI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P06682
Product Type	Transmembrane Protein
Immunogen Species	Oncorhynchus mykiss (Rainbow trout) (Salmo gairdneri)
Sequence	<p>QANGTLGRSRWLPLDPVDCVWSRWSEWTPCNSCTKIRHRSRSVEVFGQFG GKPCQGQPIG EQQRCTSDAVCEQALPSECSSIEFTCESGACIKLRLSCNGDYDCEDGSDEDC EPVRKPCG TKLYDTNEQGRTAGYGINILGMEPRINPFNNDYFNGMCNKVKNINNNEYNRLP WNVGLLN YETIAEETVSKEIYEDTYLLRELMTETKLTVSAGLNLKFTPTEKSMAKSNTTVS GGVGL DAEYDRTQMIKEVSEYTTIKNKSFMRVNGRVQLSTYRMRSRDLQVAGEFLEH VKSLPLEY EKGQYFSFLEDYGTHYTRNGKSGGEHQLVYVLNQDTIKDKKLTERRKLQDCIKV GISANFD TNIGIGGDAHIRPGHCKDTVNKNTAEKEGKALVDKVITVVRGGTLEAAVAMRT QITKEGL MDVKTYQNWARTVGDAPALLSSEPEPIQTLIPLSMPDANTRRLNMQRATQEY EAEYSVCK CKPCHNGGSLALLDGKCLCLCLPQFEGGLACQDAKADNNKNTKTPVRVFLEKG NWSCWAAW SGCSGGKRIRTRSCNTQGLSDATCRGDIVTEDYC</p>
Source	in vitro E.coli expression system
Target Names	c9
Protein Names	Recommended name: Complement component C9
Expression Region	1-574
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.