



Recombinant Human Uncharacterized aarF domain-containing protein kinase 5 (ADCK5)

Product Code	CSB-CF664749HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3MIX3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MWRPVQLCHFHSALLHRRQKPWPSPAVFFRRNVRGLPPRFSSPTPLWRKVL STAVVGAPL LLGARYVMAEAREKRRMRLVVDGMGRFGRSLKVGLQISLDYWWCTNVVLRG VEENSPGYL EVMSACHQRAADALVAGAISNGGLYVKLGQGLCSFNHLLPPEYTRTLRVLEDR ALKRGFQ EVDELFLQALPHELFQEFDYQPIAAASLAQVHRAKLHDGTSVAVKVQYIDL RDRFDG DIHTLELLLRLVEVMHPSFGFSWWLQDLKGTLAQELDFENEGRNAERCARELA HFPYVVV PRVHWDKSSKRVLTADFCAGCKVNDVEAIRSQGLAVHDIAEKLIKAFAEQIFYT GFIHSD PHPGNVLVRKGPDKAELVLLDHGLYQFLEEKDRAALCQLWRAILRDDAAMR AHAAALG VQDYLLFAEMLMQRPVRLGQLWGSHLLSREEAAYMVDMARERFEAVMAVLR ELPRPMLLV LRNINTVRAINVALGAPVDYFLMAKRAVRGWSRLAGATYRGVYGTSLLRHAK VVWEMLK FEVALRLETAMRLTALLARALVHLSLVPPAEELYQYLET
Source	in vitro E.coli expression system
Target Names	ADCK5
Protein Names	Recommended name: Uncharacterized aarF domain-containing protein kinase 5 EC= 2.7.11.-
Expression Region	1-580
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.