



# Recombinant Rat Torsin-1A-interacting protein 1 (Tor1aip1)

<b>Product Code</b>	CSB-CF728795RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q5PQX1
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Sequence</b>	<p>MAGERWRAEGLGEGWAIYVTPRAPIREGRRRLATQNGDGSAPAYETHPSR  HGRREVRFS  EPEPEVYGDFEPRAAKERSPGERRTPPEKFRSDSAKEEVRESAYNLRSRQRR  QRGPQAE  EMKTRRSTRLEQHSQQAQQQLSPATSGRGLRDAQSLSEDRGEDEPSSQPVT  SQTVSKKTV  RTPETSVMSIEDPISNLCRPPLRSPRPDASIVQHINPFEEGETEDDLESSYSDVT  IRIRSR  DSVESRDEAAVAAGHHPDSLWGLPHSRGDFTAHENQPSLLPTGCQKNPQEW  VEQAVRMRT  RMAYNNIQKSDFGNQSPSTSRQQAQVPPDESSVKIKWWLLILVAALAMGIYW  FFHTPVV  ETTAVQEFQNMKQLQSKYQSQDEKLWKRGTTFLEKHLNSSLPRPQPAILLT  AAQDAE  VLKCLSEQIADAYSSFRSVRAIRIDGAGKAAQSDLVKHEVDQELTDGFRNGQ  NAAVVHR  FESLPAGSTLIFYKYCDHENAFAFKDVALVLTVLLLEEQTLEASLGLKEIEEKVRDF  LKVKF TSSDTANSYNHMDPKLNLWSRISHLVLPVQPENALKAGSCL</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Tor1aip1
<b>Protein Names</b>	Recommended name: Torsin-1A-interacting protein 1 Alternative name(s): Lamina-associated polypeptide 1B Short name= LAP1B Lamina-associated polypeptide 1C Short name= LAP1C
<b>Expression Region</b>	1-583
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.