



Recombinant Human Suppressor of tumorigenicity 7 protein (ST7)

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-CF865116HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q9NRC1 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Homo sapiens (Human) |
| Sequence | <p>MAEAATGFLEQLKSCIVWSWTYLWTVWFFIVLFLVYILRVPLKINDNLSTVSMFL NTLTP KFYVALTGTSSLISGLILIFEWWYFRKYGTSFIEQVSVSHLRPLLGGVDNNSN NSNSSN GDSDSNRQSVSECKVWRNPLNLFARGAEYNRYTWVTGREPLTYDMNLSAQD HQTFFTCDS DHLRPADAIMQKAWRERNPQARISAAHEALEINEIRSRVEVPLIASSTIWEIKLL PKCAT AYILLAESEEATTIAEAEKLFKQALKAGDGCYRRSQQLQHHGSQYEAQHRRDTN VLVYIKR RLAMCARRLGRTREAVKMMRDLMKEFPLLSMFNIHENLLEALLELQAYADVQA VLAKYDD ISLPKSATICYTAALLKARAVSDKFSPEAASRRGLSTAEMNAVEAIHRAVEFNP HVPKYL LEMKSLILPPEHILKRGDSEAIYAFFHLAHWKRVEGALNLLHCTWEGTFRMIP YPLEKG HLFYPYPICTETADRELLPSFHEVSVYPKKELPFFILFTAGLCSFTAMLALLTHQ FPELMGVFAKAMIDIFCSAEFRDWNCKSIFMRVEDELEIPPAPQSQHFQN</p> |
| Source | in vitro E.coli expression system |
| Target Names | ST7 |
| Protein Names | Recommended name: Suppressor of tumorigenicity 7 protein Alternative name(s): Protein FAM4A1 Protein HELG |
| Expression Region | 1-585 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |