



Recombinant Danio rerio Phosphatidylinositide phosphatase SAC1-B (sacm1b)

Product Code	CSB-CF395207DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4VCH0
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	<p>MANAYERFNLHSTPEKFYIEACDDGADDVLVIDRVSTEMTLAGIKDIPPSGITRP ICGVM GTVRLVAGMYLIVITRKRKVGDLFGHTVWKAVEFDVISYKKTILHLTDIQMQDN KTFLTM INNVLNTDGFYFCTDYDLTHTQQRLSNTSPDFQEMSLLERADQRFMWNGNLL REIIAQPE LHKFAFPVIHGFIVMKPCCINGKVFEWIIISRRSCFRAGVRYYYVRGIDSEGHAAN FVETE QIVQFNNARASFVQTRGSIPFFWSQRPNLKYKPKPLISKDTNHMDGLRRHFES QVLIYGK QVILNLVNQKGSELPLEQAFKVMSSMENGFIKYIAFDHFHKECSKMRWHRLQIL VDAVSD MQEEFGYFMVSSDGKVLSEQSGTFRSNCMDCLDRTNVIQSLLARRSLQSQLQ RMGVLHVG QKIEEQADFEKIYKNAWADNANACAKQYAGTGALKTDFTRTGKRTHWGLVMD GWNSMIRY YKNNFSDGFRQDSIDLFLGNYSVDETDSLTPLHVKKDWKFLLLPVIMVVAFSM CIICLLM AGDTWTETLAYVLFWGMASALTAIVVNGREFVDAPKLVQKEKMD</p>
Source	in vitro E.coli expression system
Target Names	sacm1la
Protein Names	Recommended name: Phosphatidylinositide phosphatase SAC1-B EC= 3.1.3.- Alternative name(s): Suppressor of actin mutations 1-like protein B
Expression Region	1-586
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.