




Recombinant Varicella-zoster virus Envelope glycoprotein E (gE)

Product Code	CSB-CF872802VAQ
Abbreviation	Recombinant Varicella-zoster virus gE protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9J3M8
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Varicella-zoster virus Envelope glycoprotein E(gE)
Immunogen Species	Varicella-zoster virus (strain Oka vaccine) (HHV-3) (Human herpesvirus 3)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	SVLRYDDFHIDEDKLDTNSVYEPYYHSDHAESSWVNRGESSRKAYDHNSPYI WPRNDYDGFLENAHEHHGVYNQGRGIDSGERLMQPTQMSAQEDLGDDTGIH VIPTLNGDDRHKIVNVDQRQYGDVFKGDLNPKPQQGRLIEVSVEENHPFTLRA PIQRIYGVRYTETWSFLPSLTCTGDAAPAIQHICLKHTTCFQDVVVVDVCAENT KEDQLAEISYRFQGGKEADQPWIVVNTSTLFDLELDPPEIEPGVLKVLRTKQ YLGVIWNMRGSDGTSTYATFLVTWKGDEKTRNPTPAVTPQPRGAEFHMWN YHSHVFSVGDTFSLAMHLQYKIHEAPDLLLEWLYVPIDPTCQPMRLYSTCLYH PNAPQCLSHMNSGCTFTSPHLAQRVASTVYQNCEHADNYTAYCLGISHMEPS FGLILHDGGTTLKFVDTPESLSGLYVFVVYFNGHVEAVAYTVVSTVDHFVNAIE ERGFPTAGQPPATTKPKEITPVNPGTSPLLRYAAWTGGLAAVLLCLVIFLICT AKRMRVKAYRVDKSPYNQSMYYAGLPVDDFEDSESTDEEEFGNAIGGSHG GSSYTVYIDKTR
Research Area	Others
Source	in vitro E.coli expression system
Target Names	gE
Expression Region	31-623aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	68.4 kDa
Protein Length	Full Length of Mature Protein



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.