



Recombinant Human DBH-like monooxygenase protein 1 (MOXD1)

Product Code	CSB-CF764931HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6UVY6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	GSGRTYPHRTLLDSEGKYWLGWSQRGSQIAFRLQVRTAGYVVGFGFSPTGAM ASADIVVGG VAHGRPYLQDYFTNANRELKKDAQQDYHLEYAMENSTHTIIEFTRELHTCDIND KSITDS TVRVIWAYHHEDAGEAGPKYHDSNRGTKSLRLLNPEKTSVLSTALPYFDLVNQ DVPIPNK DTTYWCQMFKIPVFQEKKHHVIKVEPVIQRGHESLVHHILLYQCSNNFNDSVLES GHECYH PNMPDAFLT CETVIFAWAIGGEGFSYPPHVGLSLGTPLDPHYVLLLEVHYDNPT YEEGLID NSGLRLFYTMDIRKYDAGVIEAGLWVSLFHTIPPGMPEFQSEGHCTLECLEEAL EAEKPS GIHVFAVLLHAHLAGRGIRLRHFRKKGKEMKLLAYDDDFDFNFQEFQYLKEEQTI LPGDNL ITECRYNTKDRAEMTWGGLSTRSEMCLSYLLYPRINLTRCASIPDIMEQLQFI GVKEIY RPVTTWPFIIKSPKQYKNLSFMDAMNFKWTKKEGLSFNKLVLSPVNVRCCK TDNAEWS IQGMTALPPDIERPYPYKAEPLVCGTSSSSSLHRDFSINLLVCLLLL SCTLSTKSL
Source	in vitro E.coli expression system
Target Names	MOXD1
Protein Names	Recommended name: DBH-like monooxygenase protein 1 EC= 1.14.17.- Alternative name(s): Monooxygenase X
Expression Region	20-613
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.