



Recombinant Pongo abelii Torsin-1A-interacting protein 1 (TOR1AIP1)

Product Code	CSB-CF725304PYX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5R7A3
Product Type	Transmembrane Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sequence	<p>MAGEGRRAEAVREGWGVYVTPRAPIREGRGRLAPQNGGSSDAPAYRTSLSR QGRREVRFS DEPPEVYGDPEPLVDKERSPVGKRTRLEEFRSDSAKEEVRESAYYLRSRQRR QPRPQEA EMKTRRTTRLQQQHSQQPPLQPSPVMTRRGLRDSHSSEEDPSSPTDLSQTI SKKTVRSI QEAPAESEDLVISLRPPLRYPRSEATSVQQKVNFFSEEGETEDDQDSSHSSVT TVKRSR DSESGDKTTRSSSQYIESFWQSSQSQNFTAHDKQPSVLSSGYQKTPQEWA PQTARMRTR MQTSSPGKSSIIYGSFSDDDSIKSELGNQSPSTSSQQVTGQPQNASFVKRNW WWLLPLIA ALASGSFWFFSTPEVETTAVQEFQNMNQLKNKYQGQDEKLWKRSTFLEK HLNSSHPRS QPAILLLTAARDAEEALRCLSEQIADAYSSFHVSRAIRIDGTDKATQDSDTVKLE VDQEL SNGLKNGQNAAVVHRFESFPAGSTLIFYKYCDHENAFAFKDVALVLTVLLEEETL GTSLGL KEVEEKVRDFLKVKFTNSNTPNSYNHMDPKLNLWSRISHLVLPVQPENALK RGICL</p>
Source	in vitro E.coli expression system
Target Names	TOR1AIP1
Protein Names	Recommended name: Torsin-1A-interacting protein 1 Alternative name(s): Lamina-associated polypeptide 1B Short name= LAP1B
Expression Region	1-598
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.