




Recombinant Rat MICOS complex subunit Mic60 (Immt), partial

Product Code	CSB-CF665982RA
Abbreviation	Recombinant Rat Immt protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3KR86
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose
Product Type	Transmembrane Proteins
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	STSSSSGVTAGKIAGAGLLFVGGGIGGTILYAKWDSHFRESVEKTIPYSDKLFG MVLGSAPYTVPLPKKPIQSGPLKISSVSEVMTDSELPMAQTQETNGDTPASAA GDPAPEVEHEDTINTECPNTDEGTSTFVTAALAKSLEDALNQTATVTRQTITAQ NAAVQAVKAHSSTLKTAMDNSEIAGEKKSAQWRTVEGALKERRKAVDEAADA LLKAKEELEKMKTIIEDAKKREIAGATPYITAAEEKLHSMIVDLDSVVKKVQAAQ SEAKVVSQYHEL VVQARDDFRKELDSITPDITPGWKGMSSISDLAGTLSTDDL N ALIAHAHRRIDQLNRELAQQKATEKQHIELALERQKLEEKRAFD SAVAKALEHH RSEIQAEQDRKVVEVRDAMENEMRTQLRRQAAAHTDHLRDVLKVQE QELKFE FEQDLSEKLSEQELEFHRRSSEQMDNFTLDINTAYARLRGIEQAVQSHAVAEE EARKAHLWLSVEALKYSMTSSAEMPTIPLGSAVEAIRVSCSDNEFTQALTA AIPPELSTRGVYSEETLRARFYAVQKLAGRVAMIDETKNSLYQ
Research Area	Others
Source	in vitro E.coli expression system
Target Names	Immt
Protein Names	Recommended name: Mitochondrial inner membrane protein Alternative name(s): Mitofilin
Expression Region	34-609aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	66.3 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.