



# Recombinant Human Ectonucleoside triphosphate diphosphohydrolase 4 (ENTPD4)

<b>Product Code</b>	CSB-CF007693HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q9Y227
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	<p>MGRIGISCLFPASWHFSISPVGCPRILNTNLRQIMVISVLAAAVSLLYFSVVIIRN  KYGR  LTRDKKFQRYLARVTDIEATDTNPNVNYGIVVDCGSSGSRVFVYCWPRHNG  NPHDLLDI  RQMRDKNRKPVVMKIKPGISEFATSPEKVS DYISPLLNFAAEHVPRAKHKETPL  YILCTA  GMRILPESQQKAILEDLLTDIPVHFDFLFSDSHAEVISGKQEGVYAWIGINFVLG  RFEHI  EDDDEAVVEVNIPGSESSEAIVRKRTAGILDMGGVSTQIAYEVPKTVSFASSQQ  EEVAKN  LLAEFNLGCDVHQTEHVYRVYVATFLGFGGNAARQRYEDRIFANTIQKNRLLG  KQTGLTP  DMPYLDPCLPLDIKDEIQNGQTIYLRGTGDFDLCRETIQPFMKNKTNETQTSLN  GVYQPP  IHFQNSEFYGFSEFYCTEDVLRMGGDYNAAKFTKAAKDYCATKWSILRERFD  RGLYASH  ADLHRLKYQCFSAWMFEVFHRGFSFPVNYKSLKTALQVYDKEVQWTLGAIL  YRTRFLPL  RDIQQEAFRASHTHWRGVSVFYNYHYLFSGCFLVLLAILLYLLRLRRIHRRTPR  SSSAAA LWMEEGLPAQNAPGTL</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	ENTPD4
<b>Protein Names</b>	Recommended name: Ectonucleoside triphosphate diphosphohydrolase 4 Short name= NTPDase 4 EC= 3.6.1.6 Alternative name(s): Lysosomal apyrase-like protein of 70 kDa Uridine-diphosphatase Short name= UDPase
<b>Expression Region</b>	1-616
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.