



Recombinant Measles virus Hemagglutinin glycoprotein (H)

Product Code	CSB-CF357366MCQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P08362
Product Type	Transmembrane Protein
Immunogen Species	Measles virus (strain Edmonston) (MeV) (Subacute sclerose panencephalitis virus)
Sequence	<p>MSPQRDRINAFYKDNPHPKGSRIVINREHLMIDRPYVLLAVLFVMFLSLIGLLAIA GIRL HRAAIYTAEIHKSLSTNLDVTNSIEHQVKDVLTPFKIIGDEVGLRTPQRFTDLVK FISD KIKFLNPDREYDFRDLTWCINPPERIKLDYDQYCADVAAEELMNALVNSTLLET RTTNQF LAVSKGNCSGPTTIRGQFSNMSLSLLDLYLGRGYNVSSIVTMTSQGMYGGTYL VEKPNLS SKRSELSQLSMYRVFEVGVIRNPGLGAPVFHMTNYLEQPVSNDLSNCMVALG ELKLAALC HGEDSITIPYQGSGKGVSFQLVKLGVWKSPTDMQSWVPLSTDDPVIDRLYLSS HRGVIAD NQAkWAVPTTRTDDKLRMETCFQQACKGKIQCENPEWAPLKDNRIPSYGV LSVDLSLT VELKIKIASGFGPLITHGSGMDLYKSNHNNVYWLTIPPMKNLALGVINTLEWIPR FKVSP YLFNVPIKEAGEDCHAPTYLPAEVDGDVKLSSNLVILPGQDLQYVLATYDTSRV EHAVVY YVYSPSRFSYFYPFRLPIKGVPIELQVECFTWDQKLWCRHFCVLADSESGGH ITHSGME GMGVSCTVTREDGTNRR</p>
Source	in vitro E.coli expression system
Target Names	H
Protein Names	Recommended name: Hemagglutinin glycoprotein
Expression Region	1-617
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.