



Recombinant Human Uncharacterized aarF domain-containing protein kinase 2 (ADCK2)

Product Code	CSB-CF772029HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7Z695
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MVAPWRVSVRVCLSHLRFCFELRQGLSLLRPSECPRDARLCWLLLGTLPKVVS LCGDVGEG APDVLRRRVRCSCGAAGAGPAESLPRAGPLGGVFLHLRLWLRAGALLVKFFP LLLLYPLT YLAPSVSTLWLHLLLKATETSGPTYIKLGQWASTRRDLFSEAFCAQFSKLHVRV TPHPWT HTERFLRQAFGDDWGSILSFENREPVGSGCVAQVYKAYANTAFLETDSVQRL GRASCLPP FSHTGAVGGLRELFGLGNGRKPPENLADQSFLERLLLPKADLVGSNAGVSR AQVPGHQ EATNLISVAVKVLHPGLLAQVHMDLLMKIGSRVLGVLPGIKWLSLPEIVEEFK LMVQQ IDLRYEAQNLEHFQVNFRNVKAVKFPTPLRPFVTREVLVETYEESVPVSSYQQ AGIPVDL KRKIARLGINMLLKMIFVDNFVHADLHPGNILVQGANGSSSQEAQLQQADICD TLVVAV PSSLCPLRLVLLDAGIVAELQAPDLRNFRAVFMVVMGQGQQRVAELILHHARA SECRDVE GFKTEMAMLVLTQARKNTITLEKLHVSSLLSSVFKLLMTHKVKLESNFASIVFAIM VLEGL GRSLDPKLDILEAARPFLLTGPVCP
Source	in vitro E.coli expression system
Target Names	ADCK2
Protein Names	Recommended name: Uncharacterized aarF domain-containing protein kinase 2 EC= 2.7.11.-
Expression Region	1-626
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.