



Recombinant *Gossypium hirsutum* 3-hydroxy-3-methylglutaryl-coenzyme A reductase 2 (HMG2)

| | |
|--------------------------|--|
| Product Code | CSB-CF525369GHB |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | O64967 |
| Product Type | Transmembrane Protein |
| Immunogen Species | <i>Gossypium hirsutum</i> (Upland cotton) (<i>Gossypium mexicanum</i>) |
| Sequence | <p>MEASRRSAVKPVIVLKPSKRFP LVEDTPGKASDALPLPLYLTNAVFFTLFFTVV YFLLSR WREKIRASIPLHAVTFPEIVAIFAFVASLIYLLGFFGIDFVQSLIIRPSGDVWSGED YEE ENEVLLHEEDARTVPCGQALDCSVPSLPHMARNVTAQRLFDEKPVRVATEED ARKVSCGQ AVDCSLHSLPPRPPIVTSQKLFHEKTVIVTTEEDEEIIKSVVAGTLPSYSLESKLG DCKR AAAIRREALQRLTGRSLSGLPLDGFYESILGQCCEMPVGYVQIPVGIAGPLLL NGREYS VPMATTEGCLVASTNRGCKAIHLSSGATSILLKDG MTRAPVVRVRFSTAKRAAEL KFYLED P ENFDTLAVVFNRSSRFGRLQSIKCAIAGKNLYLRFT CSTGDAMGMNMVSKGV QNVLD FLQ TDFPDMDVIGISGNFCSDKKPAAVNWIEGRGKSVVCEAIIEGDVVRKVLKTSVE SLVELN MLKNLTGSAMAGALGGFNHASNIVTAIYIATGQDPAQNVESH CITMMEAVN DGKDLHI SVTMPSIEVGT VGGGTQLASQSACLNLLGVKGASKDVAGANSRMLATIVTGAV LAGELSL MSALAAGQLVKSHMKYNRSSKDMSNLSS</p> |
| Source | in vitro E.coli expression system |
| Target Names | HMG2 |
| Protein Names | Recommended name: 3-hydroxy-3-methylglutaryl-coenzyme A reductase 2 Short name= HMG-CoA reductase 2 EC= 1.1.1.34 |
| Expression Region | 1-628 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |