



Recombinant *Schizosaccharomyces pombe* Patatin-like phospholipase domain-containing protein SPAC1786.01c (SPAC1786.01c, SPAC31G5.20c)

Product Code	CSB-CF515640SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14115
Product Type	Transmembrane Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Sequence	<p>MSIPEESEINKDYTVQEDLDEFKAYTCVYKKRHDEKIEYITAQHDWNPVYEAVV PRKSKP GKDEKREGFMYPILRWPLMFTAFLCLTFVAFLYLLDRLYINCYEYFIVWRGEAR RLRKLL QEAKTYEEWKERARALDKYFGNDEWKLDPVYDYDYTLVQAVYSSLVKHRE QKDOWNALKS VLDVCVRSNFGGIDSSMLYSRTYSGTKKLVEDYVNELKVCLETVIDQRLYTAQ ERSKMFE YFSHNYGRTALCLSGGASFAIYHTGVLRLALLNQDLIPNVITGTSGGGLLAALVC TRTNEE LKQLLPELASKYQSDIGNWLDATKRYFRTGARFDEILWAKTCMYFTRGSLTF AEAYKRT GRILNISVIPSDVHSPPKLINYLTSPTVIWSAVIASCAVPGILNPIPLMTRSQSHR LIP HNFGNRFKDGSLRTDIPLSELRTQFNHFSIVSQTNPHVQVFFFSPRGTVGRP VSHRKGR GWRGGYVGS AIEQFLKYDMIKWLHVIRSLELLPRPLGTDWSSVFLQKFDGTITI WPKTKF QDFYYILSPPSVERLGYMIDAGQAATFPKLDFAARMTIEKLIKGRMMMDKPSKL GRSID GTIGTSM SRGDIEISQESASISPEDIDIVN</p>
Source	in vitro E.coli expression system
Target Names	SPAC1786.01c
Protein Names	Recommended name: Patatin-like phospholipase domain-containing protein SPAC1786.01c EC= 3.1.1.-
Expression Region	1-630
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.