



Recombinant *Arabidopsis thaliana* Ubiquitin carboxyl-terminal hydrolase 18 (UBP18)

Product Code	CSB-CF724234DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q67XW5
Product Type	Transmembrane Protein
Immunogen Species	<i>Arabidopsis thaliana</i> (Mouse-ear cress)
Sequence	<p>MHEVGFPLDLSVFTRLIATLFFLAVGVFYFLKNTAAKYFDIGAAAAGGFDRDFM AVDAED CSVCGNFSTKKCSRCKSVRYCSAECQRSDWSSGHQRNCRDYGITTLTPSAK NGLRFRASP FGDSSASSIALISERGQNKSSLKPREVLFPYEEFVEYFNWDNPELAPCGLMNC GNSCFAN VILQCLSWTRPLVAYLLEKGHKRECMRNDWCFLCEFQTHVERASQSRFPFSP MNIISRLT NIGGTLGYGRQEDAHEFMRYAIDMMQSVCLDEFGGEKIVPPRSQETTLIQYIF GLLQSQ VQCTVCNHVSDQYENMMDLIVEMHGDAGSLEECLDQFTAEEWLHGDNMYKC DRCSDYVKA CKRLTIRRAPNILTIALKRYQGGRYGKLNKRISFPETLDLNPYMSEGGDGDSDVY KLYAVI VHLDMLNASFFGHYICYIKDFCGNWYRIDDSEIESVELEDVLSQRAYMLLYSRI QARSSS SCLRSEVKDEKKTDLTDESCVKELVSSMVGAIERSSTHATIEDPVCEQSPS PSPSPS PSPSPSPSVLASECCSEVERIDTLDESNSIDDSATDHQEDVANGNKDPE VKYQAAD SWSDPTTSTPLVCTKSKPPVRDMDTKMIDAQ</p>
Source	in vitro E.coli expression system
Target Names	UBP18
Protein Names	<p>Recommended name: Ubiquitin carboxyl-terminal hydrolase 18 EC= 3.4.19.12 Alternative name(s): Deubiquitinating enzyme 18 Short name= AtUBP18 Ubiquitin thioesterase 18 Ubiquitin-specific-processing protease 18</p>
Expression Region	1-631
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.