



Recombinant *Rhizobium meliloti* RpoH suppressor (suhR)

Product Code	CSB-CF323378RKU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15715
Product Type	Transmembrane Protein
Immunogen Species	<i>Rhizobium meliloti</i> (strain 1021) (Ensifer meliloti) (Sinorhizobium meliloti)
Sequence	<p>MARTPAKYCDLVMKGGITSGIVYPNAALALARDYRFKNIGGTSAGAIAAAACAA AAVGDR RKQMKAAIAQPEERVGFEGGLAKASANLSPGFIKDLLQPAAGAGQAFRLLVTL AGNTGVL RKGVALLGSVVRIAPVETLLLLAALAGLAYAVGGQTGMIAAALPAAICAYLGGV VFAVLR IARVLRRLMGLCTGTAPDQPARRPRMVLTDWLHETLQALSGKASGQPLTFG DLWTAERY PGEPGSDRAVTLKMITTGISHQEPRSLPFESALFWFRKEFEALFPKVVVDWM VEKAGEP VTVAGEDYLLPHGADMPVLVATRMSLSFLLISAVPLHEPARRESLPGSDGE NEAEDTT SDEDEQKTVLDSTEALTTGGKKRRARPAAFRICWFSDGGISSNFPIHLFDRALP RWPTFA INLVYPETSDTGSRPEVFLPENNRQGWQRHYQPIARKSAVHELCAFVFAIVAT MQNWRDL LQSRAPGHRERIVHVSLSPQEGGLNLAMSKEVLAAVSKKGTAAGEAFARFSFE NHYWIRW RNLASALQRYTIDIAASDAYRPKIPDYEPAYALAHDATSKPPSYRFASKAEREE AARLLE KLIGEGEKWSGEGPDLTKTAPRPLPQLQIPTY</p>
Source	in vitro E.coli expression system
Target Names	suhR
Protein Names	Recommended name: RpoH suppressor
Expression Region	1-633
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.