



# Recombinant *Ajellomyces dermatitidis* Formation of crista junctions protein 1 (FCJ1)

<b>Product Code</b>	CSB-CF512908AZE
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	C5JIS0
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Blastomyces gilchristii</i> (strain SLH14081) ( <i>Blastomyces dermatitidis</i> )
<b>Sequence</b>	<p>SSTPNAAAATPELSQKATNSTSTKPPGPNDPDVRSPASPSTGSTLHPETVSKPP  QSPAVQG  QTSPGSSVQPPEHEPSPPPPRPPPAPKTGLLRKLLYLFLTGLAYAGGVWYSL  RSDNFYD  FFTEYIPYGEEAVLYLEERDFRSRFPFIARQINRRVSAPRDEGAQVMIPGRSGL  SWKVAE  EQQEASDVTKQGQHISATDANELTEETKVAEKAKEDVKSKPVAKKAAEAAEPKS  SPKVVEP  HPAKAEENTSLEAPRQPVPVAAAAIEHLGLDNEDEPVVQDLVKVFNDIITVISAD  ESASK  FSVPIAKAKEELEKIGDRIVALKNDAQESAKEEIRNAQAALDKSAAELVRHINEV  RAQDA  AEFREEFESEREKISKSYQEKVTTTELQRAHEVAEQRLRNELVEQAIELNRKFLA  DVKTLV  ENERGRLSKLAELTANVAELERLTAGWSDVIDINLRTQQQLQVAVDSVRTTLEN  SEVPRP  FIRELAAVKELASNDEVVAAAIAISISPTAYQRGIPSPAQLVDRFRRVASEVRKAS  LLPEN  AGITSHAASLVLSKVMLKKQGTPVGNDVESILTRTENLLEEGNFDEAAREMNSL  QGWAKL LSKDWLADVRRVLEVVKQALEVIETEARLRCLQVE</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	MIC60
<b>Protein Names</b>	Recommended name: Formation of crista junctions protein 1 Alternative name(s): Mitofilin
<b>Expression Region</b>	32-665
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.