



Recombinant *Saccharomyces cerevisiae* Topoisomerase I damage affected protein 7 (TDA7)

Product Code	CSB-CF466577SVH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B5VQN1
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain AWRI1631) (Baker's yeast)
Sequence	MNSNSTIGRTTLGESDTISLSFSEPSSSLNSRSTDVVFASSTLVPQQGSLTSL PPVSST ATPTYYSTSLTYDETLHTSIDVSSTSTLVSSSTDSSSSSEQDTYSSQYDPATSSY SIITPS MSIFSSTSPMSSSSSITSEWSSLTSTTPTLSSSATSLSSSWSSLSPSSLLVSS SLSLSL SSSYSDTKLFSFDSRSSIFSPSTPTVISPSYTYLSSISATSFQISTTSELSSSWFS TISS PSTTSNKDITTFPSSSRNTSTSFYSSSLSSTNDFSTISKSSKLSPSASSSTVSISTI SVPT SSSVSSSSSKVPSNRPSSSSSSDDTTSAYSSTYTFQSLQSTTSSSIPPTTQTPS TSTIST SPIPTSSQVFENTVAISSSEDSKTIYYFYTQTYDITDSSTTFVTGLPTTIAVAKSEV TSFS APSSITADMSFYQHWLDGSLDNNKNQGTSKTNTGTIVGSGVGGILICVLV VWFMLV RKRKAKRHFKEKENDSFCHEIGRRTGFPTTAQAKEASLQAQDSGSQQRNTETAS ANNPFSNE FNFKARGNPPPVPPPRNVAMNGSFQNMRSNFMQDENRFSYGSSFTYSSLG SSTQGGFST LSSNSIRLGRGLDNDISHDERNTVQNNSQGFLREII
Source	in vitro E.coli expression system
Target Names	TDA7
Protein Names	Recommended name: Topoisomerase I damage affected protein 7
Expression Region	1-636
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.